

Indications and treatment plan of transcranial magnetolectric therapy instrument 1. Scope of application:

- 1. It is used as an adjuvant therapy for functional recovery of stroke and paralysis patients.**
- 2. Can be used for adjuvant treatment of ischemic cerebrovascular disease, neurasthenia, and brain injury.**

Two. Indication:


- 1. Ischemic cerebrovascular disease ① Cerebral thrombosis and cerebral thrombotic infarction ② Lacunar infarction ③ Insufficient blood supply to the brain, brain atrophy, cerebral arteriosclerosis ④ Insufficient blood supply to the base of the vertebra.**
- 2. Brain and spinal cord injury diseases ① patients with craniocerebral injury, toxic brain injury and spinal cord injury ② patients with motor function recovery after cerebral hemorrhage ③ children with cerebral palsy and epilepsy ④ patients with brain function recovery after traumatic craniocerebral injury and intracranial tumor surgery.**
- 3. Brain functional disorders ① Parkinson's disease, depression, Alzheimer's disease, mental disorders ② Neurasthenia, insomnia, vertigo, nervous headache, brain fatigue syndrome, neurological diseases ③ Depressive neurosis, anxiety disorder, phobia.**

Ischemic cerebrovascular disease, brain injury disease, infantile cerebral palsy and limb motor dysfunction, migraine and other diseases caused by the above diseases.

Three. Treatment time and course:

Treatment time: 20 ~ 30 minutes/time, 1 ~ 2 times a day; Course of treatment: 10-20 days as a course of treatment, with an interval of 2-3 days for each course of

treatment. For chronic diseases, it is recommended to stop the treatment for 3-5 days after the end of one course of treatment and then start the next course. For patients who need long-term treatment, they need to restart the treatment at an interval of one week after 5 courses of treatment.

 The above treatment time and course of treatment are for reference only, and the doctor determines the treatment time and course of treatment according to the patient's condition.

Four. Indications and treatment plan: (see the instructions for the specific treatment site)

	Disease	Magnetic therapy function	Massage function	Lectroencephalogram bionic electric function	Body electrotherapy function	Treatment recommendations
	Infantile cerebral palsy	First weak and then strong	Optional and optional	Postauricular mastoid point	Somatic motor function The site of the obstruction	1-2 times/day A course of treatment for 10 days
	Tourette's	First weak and then strong	Optional and optional	Postauricular mastoid point	Somatic motor function The site of the obstruction	1 time/day One course of treatment is 10-20 days.
	Hyperactivity	First weak and then strong	Optional and optional	Postauricular mastoid point	Somatic motor function The site of the obstruction	2 times/day One course of treatment is 10-20 days.
	Pediatric developmental delay (speech, Intelligence, availability, physical	First weak and then strong	Optional and optional	Mastoid point behind the ear, temple	Somatic motor function The site of the	1-2 times/day A course of treatment for 10 days

	fitness, etc.)				obstruction	
	Pediatric autism	Strong gear	Optional and optional	Postauricular mastoid point	Don't choose	2 times/day 15-20 days as a course of treatment
	Brain atrophy	Weak gear	Optional and optional	Mastoid point behind the ear, temple	Somatic motor function The site of the obstruction	1 time/day 15-20 days as a course of treatment
	Insufficient cerebral blood supply	First weak and then strong	Select	Postauricular mastoid point	Somatic motor function The site of the obstruction	1-2 times/day A course of treatment for 10 days

	Cerebral thrombosis	Strong gear	Select	Postauricular mastoid point	Somatic motor function The site of the obstruction	1-2 times/day One course of treatment is 10-20 days.
	Cerebral thrombotic infarction	Strong gear	Select	Postauricular mastoid point	Somatic motor function The site of the obstruction	1-2 times/day 15-20 days as a course of treatment
	Lacunar infarction	Strong gear	Select	Postauricular mastoid point	Somatic motor function The site of the obstruction	1-2 times/day One course of treatment is 10-20 days.
	Insufficient blood supply to the base of the vertebra	Strong gear	Optional and optional	Postauricular mastoid point	Cervical vertebrae	1 time/day One course of treatment is 10-20 days.
	Cerebral arteriosclerosis	Strong gear	Select	Postauricular mastoid point	Somatic motor function The site of the obstruction	1 time/day One course of treatment is 10-20 days.
	Stroke	First weak and then strong	Optional and optional	Postauricular mastoid point	Somatic motor function The site of the obstruction	1-2 times/day One course of treatment is 10-20 days.
	Hemiplegia	First weak and then strong	Optional and optional	Postauricular mastoid point	Somatic motor function The site of the obstruction	1-2 times/day One course of treatment is 10-20 days.

	Craniocerebral injury	First weak and then strong	Don't choose	Postauricular mastoid point	Somatic motor function The site of the obstruction	1 time/day A course of treatment for 10 days
	Motor function due to brain injury Obstacle	First weak and then strong	Don't choose	Postauricular mastoid point	Somatic motor function The site of the obstruction	1-2 times/day One course of treatment is 10-20 days.
	Toxic brain injury	First weak and then strong	Optional and optional	Postauricular mastoid point	Somatic motor function The site of the obstruction	1 time/day A course of treatment for 10 days
	Spinal cord injury	First weak and then strong	Don't choose	Postauricular mastoid point	Somatic motor function The site of the obstruction	1 time/day A course of treatment for 10 days
	Epilepsy	First weak and then strong	Optional and optional	Postauricular mastoid point	Somatic motor function The site of the obstruction	1-2 times/day 15-20 days as a course of treatment
	Senile dementia	First weak and then strong	Optional and optional	Postauricular mastoid point	Somatic motor function The site of the obstruction	1-2 times/day 15-20 days as a course of treatment
	Mental disorders	First weak and then strong	Optional and optional	Postauricular mastoid point	Don't choose	1-2 times/day A course of treatment

						for 10 days
	Parkinson's disease	Strong gear	Optional and optional	Postauricular mastoid point	Somatic motor function The site of the obstruction	1 time/day A course of treatment for 10 days
	Neurasthenia (insomnia, amnesia, vertigo, nervous headache, brain fatigue syndrome, nervous disease, depressive neurosis, anxiety disorder, panic disorder, anxiety disorder, neurosis) Phobia)	First weak and then strong	Optional and optional	Mastoid point behind the ear, temple	Don't choose	For depressive neurosis, it is recommended to treat only with magnetic therapy function. 1-2 times/day One course of treatment is 10-20 days.
	Depression	First weak and then strong	Optional and optional	Optional and optional	Don't choose	In the first course of treatment, only the treatment of magnetic therapy function was carried out. 1-2 times/day A course of treatment for 10 days
	Migraine	First	Optional	Mastoid point	Optional	1-2 times/day

		weak and then strong	and optional	behind the ear, temple	cervical spine location	One course of treatment is 10-20 days.
--	--	----------------------------	-----------------	---------------------------	----------------------------	--


Indications corresponding to the 10 prescriptions for somatic electrotherapy (the default prescription for startup is 1):

- ★ **Prescription 01: a pulse current is generated in about 1 second. For mild denervation and muscle weakness. It is more suitable for the treatment of flaccid paralysis.**
- ★ **Prescription 02: For the treatment of mild denervation and myasthenia.**
- ★ **Prescription 03: For treatment of mild or moderate denervation.**
- ★ **Prescription 04: For treatment of moderate denervation.**
- ★ **Prescription 05: For the treatment of severe denervation.**
- ★ **Prescription 06: It is used for the treatment of severe or extremely severe denervation and high muscle tension (for patients with high muscle tension, low treatment intensity is appropriate at the initial stage of treatment).**
- ★ **Prescription 07: There is a strong electrical stimulation in each pulse cycle, which has the same indications as prescription 03, and is used alternately with prescription 03 in the middle and late stages of treatment.**
- ★ **Prescription 08: There is a strong electrical stimulation in each pulse**

cycle, which has the same indications as Prescription 04, and is used alternately with Prescription 04 in the middle and late stages of treatment.

★ Prescription 09: There is a strong electrical stimulation in each pulse cycle, which has the same indications as prescription 05, and is used alternately with prescription 05 in the middle and late stages of treatment.

★ Prescription 10: There is a strong electrical stimulation in each pulse cycle, which has the same indications as Prescription 06, and is used alternately with Prescription 06 in the middle and late stages of treatment.

 The indications corresponding to the above prescriptions are for reference only, and doctors can choose according to the patient's condition.