



RheinLaser

Rhein Laser Technologies Co., Ltd.

Legal Headquarters:

801, 8F, E2 Building, Future City, No.999 High-Tech Avenue,
East Lake High-Tech Development Zone, Wuhan 430206, China
(Free Trade Zone Wuhan Area)

Production Base:

C3-6/C8-2 Building, Optics Valley Union Technology City,
HuangRong district, 436032 Ezhou

Tel: +86 27 8818 9920

Email: sales@rheinlaser.com

<https://smartlasertherapy.com/>



RheinLaser

SMART ICE

Class IV Laser Therapy System

MILESTONE OF LASER THERAPY




Deliver Maximum Results with Smart Ice Therapy Laser

- 20% Faster Treatment Times
- Optimal Treatment Outcomes
- Dual/Four Wavelengths Available in One Machine
- 12-inche Large Touch Screen
- Expert Treatment Protocols
- Intelligent Cloud Management System
- Digital Thermal Imaging for Choice
- Zero Reported Side Effects



How It Works?

SMART ICE Laser Therapy System utilizes a high-powered laser to promote healing in damaged tissues. It's a drug-free, surgery-free, non-invasive pain treatment alternative.



Outcomes

Safe, effective and non-invasive

SMART ICE LASER offers scientifically-proven treatment solutions to improve patient outcomes faster.





Class IV Therapy Laser That Redefines Laser Therapy Industry

4 Wavelengths



635nm/810nm/980nm/1064nm (915nm optional)
Delivers extensive tissue stimulation and penetration to achieve optimal outcomes.



Higher Power: 30W/45W/60W




SMART ICE SERIES LASERS have the therapeutic power of 30W/45W/60W. More power, shorter treatments, faster results!

Four Treatment Modes



CW/Pulse/B/Single Mode for targeted treatment.



Features

Expert Treatment Protocols

SMART ICE LASER contains preset treatment protocols for ease of use. Even beginners can quickly operate the machine to effectively treat various conditions, eliminating the need for extensive learning curves.



Intelligent Cloud Management System

Create profiles for patients and save treatment data, standard configuration for high-end clinics.

Digital Thermal Imaging

Visualize the therapeutic effects in real-time, allowing for precise treatment adjustments.



Accessories

Attachable treatment heads

SMART ICE SERIES LASERS handpiece comes with up to 4 treatment heads, allowing you to change the width of the beam for maximum versatility.



Goggles

Our goggles are tailor-made for the wavelengths emitted by the SMART ICE SERIES LASERS and boast an optical density rating of 5+ for utmost safety, both for you and your patients.



Travel Case


The travel case is custom designed for the SMART ICE SERIES LASERS. It is waterproof, crushproof, and dustproof, ensuring unparalleled protection.






Maximize Results for Both Your Patients and Bottom Line

Wide Coverage of Indications



Enhanced Patient Satisfaction


Recovery Rate of Back Pain



*Data based on clinical trial of 120 doctors across the United States. This chart reveals the treatment outcomes that Smart Ice Laser Therapy System delivered for the common back pain. In addition, it is also highly effective in treating other diseases, such as shoulder pain, neck pain, knee pain and arthritis.

Increased Revenue Stream

How Much Revenue Are You Generating Each Year?



*Research data indicates that the clinic's annual revenue has increased by approximately 35% after the introduction of Smart Ice Laser Therapy System.

“


I wanted to extend my heartfelt thanks for introducing the Smart Ice Laser Therapy System to our clinic. The treatment effects have been outstanding, leading to significantly improved patient outcomes and soaring levels of customer satisfaction. What's even more remarkable is that, during the same period last year, our clinic's revenue has doubled. I'm honored to be able to use Smart Ice to bring relief and healing to our patients, making a real difference in their lives.

—Dr. Jason Smith from Infinite Motion Therapy Center

”



Specification



Model	SMART ICE-S60D	SMART ICE-S45D	SMART ICE-S30D	SMART ICE-S60B	SMART ICE-S30B
Max Power	60W	45W	30W	60W	30W
Wavelength	635nm+810nm+980nm+1064nm Or 635nm+810nm+915nm+980nm			810nm+980nm	
Laser Type	GaAlAs Diode Laser				
Spot Size	Φ9mm+Φ20mm+ spherical lens tube; optional Φ120mm				
Treatment Mode	CW / Pulse / Single / B Mode				
Pilot Beam	635nm, 500mW				
Cooling Method	Air-cooled				
Power Supply	AC100-240V ,50/60Hz				
Input Power	400W				
Dimension	425mmx465mmx310mm				
Weight	20kg				

About US

Industry Authority

Rhein Laser is a medical laser company focused on the R&D and application of laser technology that maximize patient outcomes and practice profitability.

No 1. Best Seller in Physical Therapy Equipment Industry on Alibaba.



CERTIFICATE



RheinLaser

 FDA-cleared technology	 Recommended by thousands	 Patent certification	 Quality assurance
 Factory price	 20 years' production experience	 Quick after-sales response	 Innovation