### Introduction

- Handy Cure s` is a portable device designed particularly for Home healthcare market
- The device is based on Quantum Medicine (QM) and combines low-level pulse laser radiation, pulse infrared radiation, visible red light and static magnetic filed
- Effective in both acute and chronic pain situations
- Speeds recovery and relief pain
- operated up to 18 hours using a rechargeable lithium battery

## How does Handy Cure s` works

### 1. Renewed Energy Production (ATP)

- Photochemical processes
- Bio-Photons → Chemical process → Stimulates renewed energy production (ATP) + enzymes
- Benefits:
- Restore impaired cellular function
- Encourage collagen and cartilages formation
- Repair tendons and ligaments
- Tissue heal and pain relief

# How does Handy Cure s` works

### 2. Improves Microcirculation

- Ischemic injury
- Effective to inflammation
- Radiances capable of:
  - Increases supply of oxygen to cells
  - Improves microcirculation
  - Remove toxic wastes
  - Tissue cured and pain relief

# How does Handy Cure s` works

### 3. Re-generation New cells

- Encourages re-generation of new cells
- Stimulate body Immune system

# Benefits of using Handy Cure s`


- Increases <u>speed</u>, <u>quality</u> and <u>strength</u> of cell tissue repair
- Increase blood supply
- Stimulates immune system
- Stimulates nerve function
- Develops collagen and muscle tissue
- Help generates new healthy cells and tissue

# Conditions treated by Handy Cure s`

- Spinal diseases comprising early osteochondrosis, scoliosis, intervertebral disk hernia, post-traumatic pain, cartilage disorders, rediculitis, ischialgia
- Joint diseases, such as arthritis, arthrosis, coxarthrosis, contractures, myositis
- Muscle and ligament strains

### **Treatment programs**

**Variable Program** 50 Hz 5 Hz Chronic **Acute** Initial **Treatment** pain pain treatment **Pulse** 1000-2500 50 Hz/Sec 5 Hz/Sec Hz/Sec frequency **Program** 5 minutes 5 minutes 5 minutes duration



### **Technical information**

#### **Average power:**

Total radiation: 60-90 mW

Laser radiation: 0.4-6.25 mW Infrared radiation: 30-90 mW Red LED's radiation: 2-10 mW

**Permanent magnet induction:** 25-45 mT

#### **Wavelength of:**

Laser radiation: 905 nm

Infrared radiation: 875 nm

Red radiation: 635 nm

Laser power: 25W

**Maximum pulsating laser power:** 25W

Class II equipment, Applied part Type B

Laser class 1 (IEC 60825-1)

Supply ratings: 13 VDC, Max. 1.3A

## Handy Cure s'

- Non-invasive treatment
- Painless
- No side effects
- Safe to use
- Easy to use
- Compact & Portable



C E 1023

